**Endo Linear Cutter Stapler  
Reloading Unit for Endo Linear Cutter Stapler**

**Endo Linear Cutter Stapler**  
**Reloading Unit for Endo Linear Cutter Stapler**  
**EN**

**⚠ BEFORE USING PRODUCT, READ THE FOLLOWING INFORMATION THOROUGHLY**  
**IMPORTANT**

This booklet is designed to assist in using this product. It is not a reference to surgical techniques. This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Do not reuse, reprocess or re-sterilize this device.

**DESCRIPTION**

The Endo Linear Cutter Stapler places two, triple-staggered rows of titanium staples and simultaneously divides the tissue between the two triple-staggered rows. The size of the staples is determined by the selection of the 2.0mm, 2.5mm, 3.5mm, 3.8mm, 4.0mm, 4.2mm or 4.4mm Reloading Unit for Endo Linear Cutter Stapler (RU). The Endo Linear Cutter Stapler will accommodate any of the reloading unit sizes that are available in the 30mm, 45mm, and 60mm stapling lines.

All Reloading units must be used with trocars of 12mm or larger.

The CEAG stapler can be used with Panther Healthcare's CADG reloading units. Each stapler can be used with reloading units up to 25 times.

**INDICATIONS**

The Endo Linear Cutter Stapler has applications in abdominal and thoracic surgery for resection, transection, and creation of anastomosis.

**CONTRAINDICATIONS**

1. The Endo Linear Cutter Staple instrument should not be used on tissue such as liver or spleen where compressibility is such that closure of the instrument would be destructive.
2. Do not use the Endo Linear Cutter Staple where adequacy of hemostasis cannot be verified visually after applications.
3. These devices are provided STERILE and are intended for use a SINGLE procedure only.

DISCARD AFTER USE. DO NOT RESTERILIZE.

**WARNINGS AND PRECAUTIONS**


1. Preoperative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected staple size. Careful consideration should be given to any pre-surgical treatment the patient may have undergone and in corresponding selection of staple size.
2. Always include the combined thickness of the tissue and of any staple line reinforcement material used when choosing the proper staple cartridge.
3. Always close the jaws of the Endo Linear Cutter Staple prior to introducing and removing the stapler from the trocar sleeve.
4. After firing, always inspect the staple line for hemostasis. Minor bleeding may be controlled by electrocautery or manual sutures.
5. Placement of tissue proximal to the tissue stops (on the RU) may result in stapler malfunction. Any tissue extending beyond the cut mark will not be transected.
6. When using the stapler more than once during a SINGLE surgical procedure, be sure to remove the empty RU and reload a new one. A safety interlock is provided that prevents an empty RU from being fired a second time. Do not attempt to override the safety interlock.
7. When positioning the stapler on the application site, ensure that no obstructions, such as clips, are incorporated into the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.
8. Endoscopic procedures should be performed only by physicians having adequate training and familiarity with endoscopic techniques. Prior to performance of any endoscopic procedures, consult the medical literature relative to techniques, complications and hazards.
9. A thorough understanding of the principles involved in laser and electro-surgical procedures is essential to avoid shock and burn hazards to both patient and operator(s), and damage to the instrument.
10. When endoscopic instruments and accessories from different manufactures are employed together in a procedure, verify compatibility, and ensure that electrical isolation or grounding is not compromised.
11. The anvil must be completely visible, (past the trocar sleeve) prior to opening the RU within the body cavity.
12. When using a staple line buttressing material, follow the instructions provided by the manufacturer of the buttress material, as performance of the stapler may be affected when using buttress materials.
13. The instrument and reloading unit are provided STERILE and are intended for use in a SINGLE procedure only.

DISCARD AFTER USE. DO NOT RESTERILIZE.

14. The device cannot be used in direct contact with the heart, the central circulatory system or the central nervous system.

① SCHEMATIC VIEW

- 1.Pin
- 2.Shaft
- 3.Non-stop articulation lever
- 4.Reverse knob
- 5.Shell
- 6.Reset button
- 7.Handle
- 8.Jaw closing trigger
- 9.Firing trigger
- 10.Reloadiing unit UNLOAD/UNLOCK BUTTON
- 11.ROTATION COLLAR
- 12.Anvil
- 13.Reloadiing unit SHAFT
- 14.Assemble interface
- 15.Knife
- 16.Cartridge holder
- 17.Cartridge
- 18. Protective cover

② LOADING

1. The RU is packaged in the open position. Do not attempt to close the RU.

**WARNING:** SELECT A RU WITH THE APPROPRIATE STAPLES SIZE FOR THE TISSUE THICKNESS. - OVERLY THICK OR THIN TISSUE MAY RESULT IN UNACCEPTABLE STAPLE FORMATION. ALWAYS INCLUDE THE COMBINED THICKNESS OF THE TISSUE AND OF ANY STAPLE LINE REINFORCEMENT MATERIAL IN USE WHEN CHOOSING THE PROPER STAPLE CARTRIDGE.

**CAUTION:** Do not attempt to remove the shipping wedge until the RU is loaded into the instrument.

2. To load the Endo Linear Cutter Stapler with appropriate RU, insert the pin located at the distal end of the instrument shaft into the RU. Ensure that the LOAD alignment indicator on the RU aligns with the LOAD alignment indicator on shaft. Push the RU in and twist clockwise 45° relative to the instrument, so that the RU will lock into place. The LOAD alignment indicator on the instrument shaft will align with the LOAD alignment indicator on the RU.


- LOAD ALIGNMENT INDICATOR (RU)
- LOAD ALIGNMENT INDICATOR (SHAFT)

3. Remove the shipping wedge from the RU prior to inserting the instrument into the trocar.

**CAUTION:** Do not close the jaw prior to removing shipping wedge.

4. Squeeze the close trigger to close the jaw and press the reset button to open the jaw )

③ UNLOADING


The knife must be in its original position and the jaw must be open when unload the reloading unit.

After firing, press down the reverse knob and squeeze the firing trigger to withdraw the knife to original position. Then press the Reset button to open the jaw.


Push down the UNLOAD/UNLOCK button and rotate the reloading unit 45 degrees counterclockwise to remove the reloading unit from the instrument.

#### ④ INSTRUCTIONS FOR USE

**NOTE:** The jaws of the RU must be closed prior to introducing the instrument into the trocar sleeve.


The instrument shaft can be rotated 360° by using the rotation collar. The CEAG instrument can control the reloading unit to rotate 60 degrees to the left and right.

Apply the Endo Linear Cutter Stapler across the tissue to be transected. Caution: Ensure that no obstructions (such as clips) are incorporated in the instrument jaws. Firing over an obstruction may result in incomplete cutting action and/or improperly formed staples.


The instrument will not cut tissue beyond the black cut mark indicated on the RU. More than one application of the Endo Linear Cutter Stapler may be necessary for tissue exceeding the length of the RU (30mm, 45mm, or 60mm).

Insert the tissue to be cut into the jaws of the device, ensuring there are no other obstructions within the jaws, as this could result in incomplete cutting or poor staple formation.

Squeeze the firing trigger until a 'click' sound be heard, indicating that the jaws are closed. Continuously squeeze the firing trigger to activate the stapler and cut the tissue until the cutting blade reaches the end of the staple cartridge base and the trigger cannot be pressed any further, signaling the completion of the firing.

**1. CAUTION: The reload unit is equipped with a safety device, and once it has been triggered, it cannot be triggered again.**

Before firing, the jaw can be opened by pressing the reset button, allowing for readjustment of the tissue clamping position.

The total number of squeezes is relative to the length of the RU (30, 45 or 60). The window on the handle displays the number of firings. For the 30mm RU, squeezing the trigger twice consecutively will display the number 2 in the window, indicating completion of the firing sequence. For the 45mm RU, squeezing the trigger three times consecutively will display the number 3 in the window, indicating completion of the firing sequence. For the 60mm

RU, squeezing the trigger four times consecutively will display the number 4 in the window, indicating completion of the firing sequence.

FAILURE TO COMPLETELY FIRE THE RU WILL RESULT IN AN INCOMPLETE CUT AND/OR INCOMPLETE STAPLE FORMATION, WHICH MAY RESULT IN POOR HEMOSTASIS.

2. After firing, press down the reverse knob ,and then squeeze the firing trigger to return the cutting knife to its initial position. Then, press the reset button to open the jaws and release the tissue. Please gently remove the instrument from the tissue. After removal, check if the suture site is hemostatic. Minor bleeding can be stopped with electrocautery or manual suturing.

**NOTE:** Do not attempt to insert or remove the instrument from the trocar sleeve if the instrument is in the articulated position.

**The Endo Linear Cutter Stapler may be reloaded and fired up to 25 times in a single procedure.**

**⑤ STAPLE SPECIFICATIONS**

T) INSTRUMENT CODES

U) RU CODES

V) STAPE LINE LENGTH

W) COLOR

X) OPEN STAPLE SIZE

Y) CLOSED STAPLE SIZE

T	U	V	W	X	Y
CEA-30 CEAG-45 CEAG-60 CEAG-30NY CEAG-45NY CEAG-60NY	CADG-30G	31mm	Brown	2.0mm	0.75mm
	CADG-30T	31mm	White	2.5mm	1.0mm
	CADG-30N	31mm	Blue	3.5mm	1.2mm
	CADG-30P	31mm	Purple	3.8mm	1.35mm
	CADG-30S	31mm	Gold	4.0mm	1.5mm
	CADG-30D	31mm	Green	4.2mm	1.8mm
	CADG-30R	31mm	Black	4.4mm	2.0mm
	CADG-45G	47mm	Brown	2.0mm	0.75mm
	CADG-45T	47mm	White	2.5mm	1.0mm
	CADG-45N	47mm	Blue	3.5mm	1.2mm
	CADG-45P	47mm	Purple	3.8mm	1.35mm
	CADG-45S	47mm	Gold	4.0mm	1.5mm
	CADG-45D	47mm	Green	4.2mm	1.8mm
	CADG-45R	47mm	Black	4.4mm	2.0mm
	CADG-60G	59mm	Brown	2.0mm	0.75mm
	CADG-60T	59mm	White	2.5mm	1.0mm
	CADG-60N	59mm	Blue	3.5mm	1.2mm
	CADG-60P	59mm	Purple	3.8mm	1.35mm
	CADG-60S	59mm	Gold	4.0mm	1.5mm
	CADG-60D	59mm	Green	4.2mm	1.8mm
CADG-60R	59mm	Black	4.4mm	2.0mm	

Benefits: Use the product, patient can have lower blood loss, shorten hospital stay time, fewer postoperative complications, faster recovery and so on.

## MR Information

	
MRI Safety Information	
A person implanted with the Disposable endoscopic linear cutting staple may be safely scanned under the following conditions. Failure to follow these conditions may result in injury.	
Device Name	Endo Linear Cutter Stapler
Static Magnetic Field Strength (B0)	1.5T or 3.0T
Maximum Spatial Field Gradient	720 gauss/cm
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	There are no Transmit Coil restrictions
Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	NA
Scan Duration	Staple rises maximum 1.5°C for 15 min of scanning as defined above
MR Image Artifact	Under T1 SE and GRE pulse sequence, Staple extended maximum 4.46mm artifact size compared to its original shape

The Implantable Staples are MR Conditional when the above conditions are met.

**STORE AT ROOM TEMPERATURE.**

**AVOID PROLONGED EXPOSURE TO ELEVATED TEMPERATURES.**

**DO NOT EXPOSE TO TEMPERATURES ABOVE 130°F (54°C).**

From the EO sterilization date, the product is valid for 5 years.

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**B.J.ZH.F.PANTHER MEDICAL EQUIPMENT CO., LTD**